Titanium® *Taq* SP DNA Polymerase GMP grade

Code No. R320A
Size: 100 rxns

Shipping on Dry ice
Store in freezer
(-25 to -15°C)

Components:

50X Titanium *Taq* SP DNA Polymerase GMP grade 100 μ l
10X Titanium *Taq* SP PCR Buffer GMP grade 600 μ l

Description:

Titanium *Taq* SP DNA Polymerase GMP grade is manufactured under rigorous standards to ensure high quality and consistency, in a facility compliant with GMP regulations and guidelines. These manufacturing practices bring an added level of confidence and consistency to this enzyme blend. It contains a 5' to 3' exonuclease-deficient *Taq* polymerase and TaqStart® Antibody, a monoclonal antibody that inhibits *Taq* and *Taq*-derived polymerases at ambient temperatures. The presence of TaqStart Antibody in the polymerase mix allows automatic hot-start PCR. Titanium *Taq* SP DNA Polymerase GMP grade is suitable to amplify targets of up to 2 kb from highly complex templates, such as mammalian genomic DNA, and up to 4 kb from a moderately abundant cDNA or a less complex genomic DNA.

Applications:

- Genotyping
- Whole genome PCR
- Multiplex PCR
- Primer extension
- Hot-start PCR
- Preparative PCR
- Rare template amplification

Primer design:

Primer design is the single largest variable in PCR applications and the single most important factor in determining the success or failure of PCR reactions. Always check your primer design before constructing or ordering primers.

Titanium *Taq* SP DNA Polymerase GMP grade can be used in a wide variety of PCR applications, and the constraints on primer design will vary from one application to the next. In general, primers should have a T_m of ~70°C to achieve optimal results in a two-step cycling program with a 68°C combined annealing/extension step. Therefore, whenever possible, primers should be at least 22 nucleotides long (25 - 30 mers are preferred) and should have a GC content of 45 - 60%. Furthermore, the 3' ends of each primer should not be complementary to each other and should have a low GC content.

General reaction mixture for PCR (50 μ l total)

Combine the following in a PCR tube on ice:

Reagent	Volume
PCR-Grade Water	40 μ l
10X Titanium <i>Taq</i> SP PCR Buffer GMP grade	5 μ l
50X dNTP Mix (10 mM ea.)	1 μ l
5' primer (10 μ M)	1 μ l
3' primer (10 μ M)	1 μ l
50X Titanium <i>Taq</i> SP DNA Polymerase GMP grade	1 μ l
DNA Template (100 ng/ μ l)	1 μ l
Total	50 μ l

Recommended cycling conditions:

Use the following guidelines when setting up your initial experiments with Titanium *Taq* SP DNA Polymerase GMP grade. These are general guidelines—the optimal cycling conditions may vary with different thermal cyclers and will depend on your particular primers, templates, and other experimental variables.

95°C	1 min	} 25 - 30 cycles
95°C	30 sec	
68°C	1 - 3 min	
68°C	3 min	
4°C		

OEM and bulk product supply:

If you are interested in Commercial Use of our products, including OEM, bulk purchasing, or custom formulation (e.g., highly concentrated reagents, custom volumes), please contact Business Development at bd_oem@takarabio.com to discuss your needs, or visit http://info.takarabio.com/OEM_inquiry

Titanium and TaqStart are registered trademarks of Takara Bio USA, Inc.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from TAKARA BIO INC.

If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at www.takara-bio.com.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

Titanium® Taq SP DNA Polymerase GMP grade

Code No. R320A
Size : 100 rxns

Shipping on Dry ice
Store in freezer
(-25 to -15°C)

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50×Titanium Taq SP DNA Polymerase GMP grade 100 μl
10×Titanium Taq SP PCR Buffer GMP grade 600 μl

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