# Anti-Human Fibronectin, Monoclonal (Clone FN30-8)

Code No. M010

Size : 0.4 mg Mouse Ig

Subclass : IgG1

\* 2 years from date of receipt under proper storage conditions.

## Source :

Monoclonal antibody was obtained by fusing the mouse myeloma cell-line P3U1 with spleen cells of BALB/c mouse after immunization with human plasma fibronectin. The monoclonal antibody was harvested from ascitic fluid.

## Purification :

Antibody was purified by column chromatography, dissolved in 10 mM PBS, pH 7.4, containing 1.0% bovine serum albumin, and then lyophilized. The lyophilized antibody does not contain preservative.

Form : Lyophilized

## Reconstitution :

Dissolve the lyophilized antibody in 200  $\mu$ l of distilled water (final concentration 2.0 mg/ml). This solution can be used as a stock solution. If dilution is necessary for your application, dilute the stock solution with the following Dilution solution just prior to use. When the entire amount of antibody is to be used over a short time period, it may be dissolved directly in 500  $\mu$ l or more of the Dilution solution.

Note (1) : Be sure to store the antibody at a minimum concentration of 2.0 mg/ml. A lower antibody concentration may result in decreased stability.

Note (2) : Reconstituted antibody solution should contain 0.1% sodium azide as a preservative when stored at 4°C.

**Dilution solution :** 10 mM PBS (pH 7.4)  
1.0% BSA  
(0.1% NaN<sub>3</sub>) \*

\* When stored at 4°C, 0.1% sodium azide should be added as a preservative.

## Working concentration:

- Western blotting :  
1 - 10  $\mu$ g/ml (different in sensitivity of detection system)
- Immunohistochemical staining :  
10  $\mu$ g/ml
- Inhibition Assay :  
final 0.1 - 1  $\mu$ g/ml (on 10  $\mu$ g/ml Fibronectin coating plate)

**Storage :** 4°C

This product does not contain preservative.

The stock solution (2.0 mg/ml) should be stored in aliquots at -20°C for 1 year, or should be stored at 4°C for 6 months after adding 0.1% sodium azide. Avoid repeated freeze-thaw cycles. Diluted antibody should not be stored.

## Application:

- Analysis of fibronectin domain structure and cell adhesion mechanism<sup>1,3)</sup>
- Adhesion blockade studies against fibronectin
- Western blot analysis under reducing or non-reducing conditions<sup>1,2)</sup>
- Immunohistochemical detection of fibronectin on frozen or paraffin embedded tissue sections \*  
\* Antigen retrieval: Proteinase K treatment

## References:

- 1) Katayama M, Hino F, Odate Y, Goto S, Kimizuka F, Kato I, Titani K, and Sekiguchi K. *Exp Cell Res.* (1989) **185**: 229-236.
- 2) Katayama M, Hino F, Kamihagi K, Sekiguchi K, Titani K, and Kato I. *Clin Chem.* (1991) **37**: 466-471.
- 3) Matsumoto Y, Sakai I, Makabe T, Yoneda J, Murata J, Kimizuka F, Ishizaki Y, Kato I, and Azuma I. *Jpn J Cancer Res.* (1991) **82**: 1130-1138.

## Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc. If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at [www.takara-bio.com](http://www.takara-bio.com). Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

# Anti-Human Fibronectin, Monoclonal (Clone FN30-8)

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
## Specificity:

Cat.#	Clone No.	Domain specificity	Inhibition of Cell Adhesion fibronectin	Cross with bovine fibronectin	Cross with porcine fibronectin	Cross with rabbit fibronectin	Cross with rat fibronectin
M001	FN8-12	Second Fibrin	None	None	None	None	None
M002	FN12-8	Cell	Yes	Yes	None	None	None
M009	FN9-1	First Fibrin-Heparin	None	Yes	Yes	Yes	None
M010	FN30-8	Cell	Yes	None	None	None	None
M013	FN1-1	C-terminal S-S bond	None	Yes	Yes	Yes	Yes (weak)
M014	FN21-1	Junction of Collagen-Heparin	None	None	Yes	Yes	None
M114	FNC4-4	Collagen	None	Yes	Yes	Yes	None
M115	FNH3-8	Type III (H12) in third Heparin	None	None	None	Yes	None

Illustration of limited digestion fragments by thermolysin and the domains recognized by each monoclonal antibody.

- 1: fibrin-heparin binding domain
- 2: collagen binding domain
- 3: heparin binding domain
- 4: cell binding domain
- 5: heparin binding domain
- 6: fibrin binding domain

SAS: Synergistic adhesion site  
 RGDS: Arg-Gly-Asp-Ser  
 CS: Connecting segment